



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

MSN:Tauto_shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0464	Structure www.ChemDrug.com
PRODUCT	Fraxinol/白蜡树精	
STANDARD TYPE	Quantitative Reference Material (含量)	
LOT NUMBER	15100922	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2015-10-09	
DATE OF REPORT	2015-10-09	

CHEMICAL NAME 2H-1-Benzopyran-2-one, 6-hydroxy-5,7-dimethoxy-

OTHER NAME

CHEMICAL FORMULA C₁₁H₁₀O₅

MOLECULAR WEIGHT(MW) 222.2

CAS NUMBER 486-28-2

CHEMICAL FAMILY 香豆素类

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK600—2008	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK600—2008	NA	ND
Water	Q/TFFK600—2008	≤2%	Conforms
Appearance		Pale White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst

Checker

Manager

Product Analysis Report

Product Name:

Fraxinol

Lot Number:

15100922

Detector:

UV

Detect Time:

2015-10-09

Analyst:

shuzhilin

Chrom-Column:

C18

Mobile Phase:

ACN: H2O (30:70)

Wavelength:


254nm

Inject Vol:

20ul

Flow:

1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	15.589	362006.438	8750307.000	100.0000