



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0350	Structure 
PRODUCT	Indolol / 吲哚醇 / 6-羟基吲哚	
STANDARD TYPE	Quantitative Reference Standard/ 含量测定	
LOT NUMBER	14082932	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2014-08-29	
DATE OF REPORT	2014-08-29	

CHEMICAL NAME	1H-indol-6-ol
OTHER NAME	6-Hydroxyindole
CHEMICAL FORMULA	C ₈ H ₇ NO
MOLECULAR WEIGHT(MW)	133.15
CAS NUMBER	2380-86-1
CHEMICAL FAMILY	Alkaloids生物碱类
FROM	

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK355—2007	NLT ≥96.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK355—2007	NA	ND
Water	Q/TFFK355—2007	≤2%	Conforms
Appearance		Pale Red Powder	Conforms

STORAGE CONDITIONS

STORAGE	Refer to the label
SHELF LIFE	2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Indanol Chrom-Column: Arcus C18
Lot Number: 14082932 Mobile Phase: ACN:H2O (Grads 20%-40%)
Detector: UV Wavelength: 231nm
Detect Time: 2014-08-29 Inject Vol: 10ul
Analyst: xiayuexia Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	6.802	258.511	16674.340	0.1581
2	8.908	302336.031	10245122.000	97.1167
3	11.520	13712.114	287493.469	2.7252