



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com

export2@tautobiotech.com

PRODUCT ID	E-0668	<p align="center">Structure</p>
PRODUCT	Cornuside/山茱萸新苷	
STANDARD TYPE	Qualitative Reference Material/ 含量测定	
LOT NUMBER	14052922	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2014-05-29	
DATE OF REPORT	2014-05-29	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₂₄H₃₀O₁₄

MOLECULAR WEIGHT(MV)

542.49

CAS NUMBER

131189-57-6

CHEMICAL FAMILY

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK522—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK522—2007	NA	ND
Water	Q/TFFK522—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Cornuside **Chrom-Column:** C18
Lot Number: 14052922 **Mobile Phase:** ACN:H₂O=22:78
Detector: UV **Wavelength:** 204nm
Detect Time: 2014-5-29 **Inject Vol:** 20ul
Analyst: shuzhilin **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	4.483	670.997	7420.049	0.1109
2	7.372	9.378	8315.306	0.1243
3	8.793	694.234	29049.602	0.4341
4	11.347	269885.313	6643440.500	99.2780
5	14.915	154.757	3532.351	0.0528