



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

MSN:Tauto_shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0755	Structure 
PRODUCT	Dehydroandrographolide Succinate / 脱水穿心莲内酯琥珀酸半酯	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	15081823	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2015-08-18	
DATE OF REPORT	2015-08-18	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₂₈H₃₆O₁₀

MOLECULAR WEIGHT(MW)

532.58

CAS NUMBER

786593-06-4

CHEMICAL FAMILY

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK597—2008	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK597—2008	NA	ND
Water	Q/TFFK597—2008	≤2%	Conforms
Appearance		White powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Dehydroandrographolide succinate Chrom-Column: C18
Lot Number: 15081823 Mobile Phase: ACN:0.1%H3PO4(grade)
Detector: UV Wavelength: 254nm
Detect Time: 2015-08-18 Inject Vol: 20ul
Analyst: shuzhilin Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	5.009	522.636	4250.300	0.0756
2	6.652	653510.688	5589757.000	99.4193
3	8.536	824.710	7562.451	0.1345
4	8.711	365.387	3009.377	0.0535
5	8.969	952.387	9322.971	0.1658
6	15.906	859.420	8506.500	0.1513

	grads	
Time/min	0	30
ACN	50	95
0.1%H3PO4	50	5