



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0782	Structure 
PRODUCT	Alisol C monoacetate/泽泻醇 C-23-醋酸酯	
STANDARD TYPE	Qualitative Reference Material/ 含量测定	
LOT NUMBER	15071321	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2015-07-13	
DATE OF REPORT	2015-07-13	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₃₂H₄₈O₆

MOLECULAR WEIGHT(MW)

528.73

CAS NUMBER

26575-93-9

CHEMICAL FAMILY

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK460—2007	NLT ≥96.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK460—2007	NA	ND
Water	Q/TFFK460—2007	≤2%	Conforms
Appearance		White Crystal	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst


Checker

Manager

Product Analysis Report

Product Name: Alisol C monoacetate
Lot Number: 15071321
Temperature: 25°C
Detect Time: 2015-07-13
Analyst: shuzhilin

Chrom-Column: C18
Mobile Phase: ACN:H2O(55:45)
Wavelength: 208nm
Inject Vol: 20ul
Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	10.887	181.154	3275.882	0.2003
2	13.622	74574.375	1574365.375	96.2837
3	14.422	527.804	10147.278	0.6206
4	17.574	292.350	9110.700	0.5572
5	21.709	392.768	12299.700	0.7522
6	24.661	387.968	14588.423	0.8922
7	29.663	194.432	11343.600	0.6937