



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0723	Structure 
PRODUCT	Cimigenol-3-O-beta-D-xylpyranoside/升麻醇-3-O-beta-D-吡喃木糖苷	
STANDARD TYPE	Quantitative Reference Standard/含量测定	
LOT NUMBER	15122922	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2015-12-29	
DATE OF REPORT	2015-12-29	

CHEMICAL NAME

OTHER NAME

Cimigenol-3-O-beta-D-xylpyranoside

CHEMICAL FORMULA

C35H56O9

MOLECULAR WEIGHT(MW)

620.81

CAS NUMBER

27994-11-2

CHEMICAL FAMILY

Glycosides 苷类

FROM

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK738-2010	NLT $\geq$ 98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK738-2010	NA	ND
Water	Q/TFFK738-2010	$\leq$ 2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst

Checker

Manager

Product Analysis Report

Product Name: Cimigenol-3-O- β -D-xylpyranoside

Lot Number: 15122922

Detect Time: 2015/12/29

Inject Vol: 20 μ L

Flow: 1.0ml/min

Analyst: zhufangfang

Chrom-Column: C18


Mobile Phase: ACN:H2O=58:42

Detector: ELSD

Tube Temp: 100 $^{\circ}$ C

Gas Flow: 2.5L/min

Impactor: Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	6.257	237037.141	2335242.000	100.0000