


PRODUCT ID	E-0691	<p align="center">Structure</p> 
PRODUCT	Specnuezhenide/特女贞苷	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	16070422	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2016-07-04	
DATE OF REPORT	2016-07-04	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₃₁H₄₂O₁₇

MOLECULAR WEIGHT(MV)

686.62

CAS NUMBER

39011-92-2

CHEMICAL FAMILY

FROM

木犀科植物女贞Ligustrum lucidum Ait. 的干燥成熟果实

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK538—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK538—2007	NA	ND
Water	Q/TFFK538—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Specnuezhenide **Chrom-Colun** C18
Lot Number: 16070422 **Mobile Phase:** MeOH:H2O=40:60
Detector: UV **Wavelength:** 224nm
Detect Time: 2016-07-04 **Inject Vol:** 20ul
Analyst: shuzhilin **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	8.252	242.054	3397.200	0.0759
2	10.203	1096.046	23218.855	0.5188
3	10.587	862.186	13186.614	0.2947
4	11.671	241186.281	4428265.000	98.9511
5	14.022	354.629	7137.800	0.1595