



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0872	Structure 
PRODUCT	Deapio platycodin D/去芹菜糖桔梗皂苷D/去芹糖桔梗皂苷D	
STANDARD TYPE	Quantitative Reference Standard/含量测定	
LOT NUMBER	16122922	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2016-12-29	
DATE OF REPORT	2016-12-29	

CHEMICAL NAME

OTHER NAME

Deapio platycodin D

CHEMICAL FORMULA

C₅₂H₈₄O₂₄

MOLECULAR WEIGHT(MW)

1093.22

CAS NUMBER

78763-58-3

CHEMICAL FAMILY

Triterpenoids三萜类

FROM

桔梗科植物桔梗 (Platycodon grandiflorus)

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK502—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK502—2007	NA	ND
Water	Q/TFFK502—2007	≤2%	Conforms
Appearance		White powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Deapio platycodin D Chrom-Column: C18
Lot Number: 16122922 Mobile Phase: 0.1%HAC: ACN (72: 28)
Detect Time: 2016-12-29 Detector: ELSD
Inject Vol: 20 µL Tube Temp: 105°C
Flow: 1.0ml/min Gas Flow: 2.9 L/min
Analyst: zhanghong Impactor: Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	8.340	647227.500	10880126.000	100.0000