



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

MSN:Tauto_shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0855	Structure 
PRODUCT	Ginsenoside RH4 /人參皂昔RH4	
STANDARD TYPE	Qualitative Analysis/鉴别用	
LOT NUMBER	16102824	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2016-10-28	
DATE OF REPORT	2016-10-28	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₃₆H₆₀O₈

MOLECULAR WEIGHT(MW)

620.43

CAS NUMBER

174721-08-5

CHEMICAL FAMILY

Triterpenoids三萜类

FROM

人參Panax ginseng C.A.Mey

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK714—2010	NLT ≥94.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK714—2010	NA	ND
Water	Q/TFFK714—2010	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Ginsenoside RH4 Chrom-Column: C18
Lot Number: 16102824 Mobile Phase: ACN: 1%H₃PO₄=45:55
Temperature: 25°C Wavelength: 203nm
Detect Time: 2016-10-28 Inject Vol: 20μL
Analyst: shuzhilin Flow: 1ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	4.227	2400.061	22384.004	2.6443
2	4.358	884.500	5617.164	0.6636
3	5.575	387.338	3755.500	0.4437
4	14.628	45866.840	802362.188	94.7869
5	15.200	396.708	6445.557	0.7614
6	16.072	253.249	5926.530	0.7001