



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0392	Structure 
PRODUCT	Jujuboside B/酸枣仁皂苷B	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	17021722	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2017-02-17	
DATE OF REPORT	2017-02-17	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₅₂H₈₄O₂₁

MOLECULAR WEIGHT(MW)

1045.21076

CAS NUMBER

55466-05-2

CHEMICAL FAMILY

Triterpenoids三萜类

FROM

鼠李科植物酸枣Ziziphus psinosa Hu的干燥成熟种子

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK234—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK234—2007	NA	ND
Water	Q/TFFK234—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Jujuboside B Chrom-Column: C18
Lot Number: 17021722 Mobile Phase: MeOH: H2O (65: 35)
Detect Time: 2017-2-17 Detector: ELSD
Inject Vol: 20 μ L Tube Temp: 105°C
Flow: 1.0ml/min Gas Flow: 2.5 L/min
Analyst: zhanghong Impactor: Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	15.857	350264.094	20350814.000	100.0000