



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0016	Structure 
PRODUCT	Ginkgolic Acid (C13:0) /银杏酸(C13:0)	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	17021423	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2017-02-14	
DATE OF REPORT	2017-02-14	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₂₀H₃₂O₃

MOLECULAR WEIGHT(M

320.4698

CAS NUMBER

20261-38-5

CHEMICAL FAMILY

Phenolic Acid酚酸类

FROM

银杏科银杏Ginkgo biloba外种皮

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK33—2005	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK33—2005	NA	ND
Water	Q/TFFK33—2005	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE
SHELF LIFE

Refer to the label
2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Ginkgolic Acids13:0 Chrom-Column: C18
Lot Number: 17021423 Mobile Phase: MeOH:0.1%H3PO4=90: 10
Temperature: 45℃ Wavelength: 310nm
Detect Time: 2017-01-14 Inject Vol: 20μL
Analyst: zhanghong Flow: 1ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	18.457	85917.844	3737539.750	98.8651
2	20.165	916.860	42906.191	1.1350