



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0015/E-0954	Structure
PRODUCT	Ginkgolic Acids (C13:0,C15:1,C17:2,C15:0,C17:1) / 总 银 杏 酸 (C13:0,C15:1,C17:2,C15:0,C17:1)	
STANDARD TYPE	Qualitative Analysis/鉴别用	
LOT NUMBER	17011106	
ASSAY METHOD	HPLC	
DATE OF SAMPLE	2017-1-11	
DATE OF REPORT	2017-1-11	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

MOLECULAR WEIGHT(MW)

CAS NUMBER

CHEMICAL FAMILY

Phenolic Acid酚酸类

FROM

银杏科银杏Ginkgo biloba外种皮

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK35—2005	NLT ≥99.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK35—2005	NA	ND
Water	Q/TFFK35—2005	≤2%	Conforms
Appearance		Pale Yellow Powder	Conforms

STORAGE CONDITIONS

STORAGE
SHELF LIFE

Refer to the label
2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Ginkgolic Acids Chrom-Column: C18
Lot Number: 17011106 Mobile Phase: MeOH:0.1%H3PO4=90:10
Temperature: 45°C Wavelength: 310nm
Detect Time: 2017-01-11 Inject Vol: 20μL
Analyst: zhanghong Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	12.440	293.356	7366.900	0.0397
2	17.957	136.500	4700.000	0.0253
3	21.890	53504.574	2345785.250	12.6270
4	23.657	242195.391	11997601.000	64.5813
5	27.723	8153.926	441639.625	2.3773
6	30.023	122.333	8452.477	0.0455
7	33.090	52.462	3873.600	0.0209
8	36.157	4910.244	343122.938	1.8470
9	38.890	30289.725	3424979.250	18.4361