



**同田生物**  
**TAUTO BIOTECH**

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

MSN:Tauto\_shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

<b>PRODUCT ID</b>	E-0484	<b>Structure</b> 
<b>PRODUCT</b>	Ciwujianoside B/刺五加皂苷B	
<b>STANDARD TYPE</b>	Quantitative Reference Material (含量测定)	
<b>LOT NUMBER</b>	17070621	
<b>ASSAY METHOD</b>	HPLC,MS,NMR	
<b>DATE OF SAMPLE</b>	2017-07-06	
<b>DATE OF REPORT</b>	2017-07-06	

**OTHER NAME**

**CHEMICAL FORMULA**

C<sub>58</sub>H<sub>92</sub>O<sub>25</sub>

**MOLECULAR WEIGHT(MW)**

1189.35

**CAS NUMBER**

114902-16-8

**CHEMICAL FAMILY**

Triterpenoids三萜类

**FROM**

刺五加Acanthopanax senticosus (Rupr. Maxim.) Harms

**ANALYTICAL RESULT**

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK136—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK136—2007	NA	ND
Water	Q/TFFK136—2007	≤2%	Conforms
Appearance		White Powder	Conforms

**STORAGE CONDITIONS**

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst

Checker

Manager

# Product Analysis Report


Product Name: Ciwujianoside B Chrom-Column: C18

Lot Number: 17070621 Mobile Phase: MeOH: H2O=25:75

Detector: UV Wavelength: 270nm

Detect Time: 2017-07-06 Inject Vol: 20ul

Analyst: shuzhilin Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	9.311	191.636	4784.700	0.0653
2	11.963	259907.000	7309944.500	99.7474
3	13.330	170.800	4731.900	0.0646
4	16.065	280.538	8998.839	0.1228