



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0130	Structure 
PRODUCT	Astragaloside II /黄芪皂苷II	
STANDARD TYPE	Quantitative Reference Standard/ 含量测定	
LOT NUMBER	17051731	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2017-05-17	
DATE OF REPORT	2017-05-17	

OTHER NAME

CHEMICAL FORMULA

C₄₃H₇₀O₁₅

MOLECULAR WEIGHT(MW)

827.011

CAS NUMBER

84676-89-1

CHEMICAL FAMILY

Saponins 皂苷类

FROM

Astragalus membranaceus Bge. Var *mongholicus* (Bge) Hsiao

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK134—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK134—2007	NA	ND
Water	Q/TFFK134—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst

Checker

Manager

Product Analysis Report

Product Name: Astragaloside II **Chrom-Column:** Arcus C18


Lot Number: 17051731 **Mobile Phase:** ACN: H2O=35:65

Detect Time: 2017-5-17 **Detector:** ELSD

Inject Vol: 20μL **Tube Temp:** 105℃

Flow: 1.0ml/min **Gas Flow:** 2.5L/min

Analyst: zhufangfang **Impactor:** Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	10.657	94292.219	1389158.625	100.0000