



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com

export2@tautobiotech.com

PRODUCT ID	E-0696	Structure 
PRODUCT	Macranthoidin B / 灰毡毛忍冬皂 昔乙	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	17042931	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2017-04-29	
DATE OF REPORT	2017-04-29	

CHEMICAL NAME

OTHER NAME

CHEMICAL FORMULA

C₆₅H₁₀₆O₃₂

MOLECULAR WEIGHT(MV)

1399.52

CAS NUMBER

136849-88-2

CHEMICAL FAMILY

FROM

忍冬科忍冬属植物灰毡毛忍冬 *Lonicera macranthoides*
Hand.-Mazz. 的干燥花蕾

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK760—2010	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK760—2010	NA	ND
Water	Q/TFFK760—2010	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Macranthoidin B **Chrom-Column:** Arcus C18
Lot Number: 17042931 **Mobile Phase:** MeOH:H2O (70:30)
Detect Time: 2017-4-29 **Detector:** ELSD
Inject Vol: 20.0 μ L **Tube Temp:** 100 $^{\circ}$ C
Flow: 1.0ml/min **Gas Flow:** 2.1 L/min
Analyst: zhufangfang **Impactor:** Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	6.690	873951.813	18755880.000	100.0000