



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

<b>PRODUCT ID</b>	E-3019	<b>Structure</b>  
<b>PRODUCT</b>	Endothal/ 去甲斑蝥酸	
<b>STANDARD TYPE</b>	Quantitative Reference Standard/ 含量测定	
<b>LOT NUMBER</b>	17081124	
<b>ASSAY METHOD</b>	HPLC,MS,NMR	
<b>DATE OF SAMPLE</b>	2017-08-11	
<b>DATE OF REPORT</b>	2017-08-11	

<b>OTHER NAME</b>	草多索
<b>CHEMICAL FORMULA</b>	C <sub>8</sub> H <sub>10</sub> O <sub>5</sub>
<b>MOLECULAR WEIGHT(MW)</b>	186.16
<b>CAS NUMBER</b>	145-73-3
<b>CHEMICAL FAMILY</b>	Monoterpenoids单萜类
<b>FROM</b>	黄黑小斑蝥Mylabrischichorri L.

#### ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK107—2007	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK107—2007	NA	ND
Water	Q/TFFK107—2007	≤2%	Conforms
Appearance		White Powder	Conforms

#### STORAGE CONDITIONS

STORAGE	Refer to the label
SHELF LIFE	2 Years


Analyst

Checker

Manager

# Product Analysis Report

**Product Name:** Endothal      **Chrom-Column:** C18  
**Lot Number:** 17081124      **Mobile Phase:** ACN:(0.1%TEA+0.1%  
H3PO4)=8:92  
**Detector:** UV      **Wavelength:** 210nm  
**Detect Time:** 2017-08-11      **Inject Vol:** 20ul  
**Analyst:** shuzhilin      **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	9.252	63573.082	831687.563	98.5478
2	17.785	511.643	12255.448	1.4522