



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0183	Structure 
PRODUCT	Ginsenoside-Rb2 /人参皂苷Rb2	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	17121221	
ASSAY METHOD	HPLC,NMR	
DATE OF SAMPLE	2017-12-12	
DATE OF REPORT	2017-12-12	

OTHER NAME

CHEMICAL FORMULA

C₅₃H₉₀O₂₂

MOLECULAR WEIGHT(MW)

1079.2685

CAS NUMBER

11021-13-9

CHEMICAL FAMILY

Triterpenes 三萜类

FROM

Panax ginseng C.A.Mey

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK76—2006	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK76—2006	NA	ND
Water	Q/TFFK76—2006	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst

Checker

Manager

Product Analysis Report


Product Name: Ginsenoside-Rb2 **Chrom-Column:** C18

Lot Number: 17121221 **Mobile Phase:** ACN:H2O=32:68

Temperature: 25℃ **Wavelength:** 203nm

Detect Time: 2017-12-12 **Inject Vol:** 20μL

Analyst: shuzhilin **Flow:** 1ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	13.955	595.449	12365.159	0.3935
2	23.525	74326.672	3096461.750	98.5472
3	25.303	310.656	12027.087	0.3828
4	27.217	438.152	18116.316	0.5766
5	28.960	68.487	3141.455	0.1000