


Manager

Product Analysis Report

Product Name: Kaempferol-3-O-rutinoside **Chrom-Column:** C18
Lot Number: 17071421 **Mobile Phase:** ACN:H2O=18:82
Detector: UV **Wavelength:** 254nm
Detect Time: 2017-07-14 **Inject Vol:** 20ul
Analyst: shuzhilin **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	9.795	323.424	5835.200	0.0897
2	13.947	209.347	5150.100	0.0791
3	16.865	278.500	7069.600	0.1086
4	19.000	164.697	4425.300	0.0680
5	23.268	177670.984	6474616.500	99.5049
6	28.771	215.400	9732.100	0.1496