


<b>PRODUCT ID</b>	E-0036	<b>Structure</b> 
<b>PRODUCT</b>	Huperzine A /石杉碱甲	
<b>STANDARD TYPE</b>	Quantitative Reference Standard/ 含量测定	
<b>LOT NUMBER</b>	17062816	
<b>ASSAY METHOD</b>	HPLC,MS,NMR	
<b>DATE OF SAMPLE</b>	2017-06-28	
<b>DATE OF REPORT</b>	2017-06-28	

<b>OTHER NAME</b>	Fordine, Selagine
<b>CHEMICAL FORMULA</b>	C <sub>15</sub> H <sub>18</sub> N <sub>2</sub> O
<b>MOLECULAR WEIGHT(MW)</b>	242.31622
<b>CAS NUMBER</b>	102518-79-6
<b>CHEMICAL FAMILY</b>	Alkaloids生物碱类
<b>FROM</b>	石杉科植物千层塔Huperzina serrata(Thumb.)Trev

#### ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK3—2006	NLT ≥98.0%	Conforms
NMR		Conforms	Conforms
MS		Conforms	Conforms
<b>Residual Solvent</b>	Q/TFFK3—2006	NA	ND
<b>Water</b>	Q/TFFK3—2006	≤2%	Conforms
<b>Appearance</b>		White Powder	Conforms

#### STORAGE CONDITIONS

STORAGE	Refer to the label
SHELF LIFE	2 Years

Analyst

Checker

Manager

