

## Manager

# Product Analysis Report

Product Name: 20(S)-Ginsenoside F<sub>2</sub> Chrom-Column C18

Lot Number: 18030925

Mobile Phase: ACN:0.1%H<sub>3</sub>PO<sub>4</sub>(40:60)

Detector: UV


Wavelength: 203nm

Detect Time: 2018-03-09

Inject Vol: 20ul

Sample dissolution: MeOH

Flow: 1.0ml/min



| No. | Ret.Time<br>min | Height<br>mAU | Area<br>mAU*min | Rel.Area<br>% |
|-----|-----------------|---------------|-----------------|---------------|
| 1   | 9.435           | 93.759        | 1229.900        | 0.0117        |
| 2   | 9.998           | 171.845       | 2504.004        | 0.0238        |
| 3   | 10.605          | 170.475       | 4645.378        | 0.0441        |
| 4   | 11.930          | 4824.602      | 115748.742      | 1.0996        |
| 5   | 13.597          | 108.357       | 2599.595        | 0.0247        |
| 6   | 15.857          | 273893.125    | 10334559.000    | 98.1792       |
| 7   | 17.868          | 51.134        | 1903.156        | 0.0181        |
| 8   | 18.798          | 39.438        | 1041.500        | 0.0099        |
| 9   | 26.577          | 1298.865      | 61988.406       | 0.5889        |