

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线: 400-820-3059

86-21-51320588-8008/8035/8026/8025

传真: 86-21-51320502

QQ:979263790

MSN:Tauto shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0664	Structure
PRODUCT	Wogonoside/汉黄芩苷	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	18010635	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-01-06	
DATE OF REPORT	2018-01-06	

OTHER NAME

Oroxindin; Glychionide B

CHEMICAL FORMULA

C22H20O11

MOLECULAR WEIGHT(MW)

460 4

CAS NUMBER

51059-44-0

CHEMICAL FAMILY

Flavonoids 黄酮类

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK542—2007	NLT $\geq$ 98.0%	99.5%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK542—2007	NA	ND
Water	Q/TFFK542—2007	$\leq$ 2%	Conforms
Appearance		pale yellow powder	Conforms

STORAGE CONDITIONS

STORAGE

Refer to the label

SHELF LIFE

2 Years

Analyst

Checker

Manager

Product Analysis Report

Product Name: Wogonoside

Chrom-Column: Kromasil C18

Lot Number: 18010635

Mobile Phase: MeOH:0.1%H3PO4=6
0:40

Detector: UV


Wavelength: 275nm

Detect Time: 2018-01-06

Inject Vol: 20ul

Analyst: zhanghong

Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	7.707	229.983	3433.295	0.0779
2	8.398	243697.828	4386922.000	99.5635
3	10.265	718.846	15799.800	0.3586