



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-3192	<p style="text-align: center;">Structure</p>
PRODUCT	Asiaticoside B/积雪草苷B	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	18041322	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-04-13	
DATE OF REPORT	2018-04-13	

OTHER NAME NA

CHEMICAL FORMULA C₄₈H₇₈O₂₀

MOLECULAR WEIGHT(MW) 975.12

CAS NUMBER 125265-68-1

CHEMICAL FAMILY Triterpenoids

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK20—2004	NLT ≥98.0%	100.0%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK20—2004	NA	ND
Water	Q/TFFK20—2004	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE Refer to the label

SHELF LIFE 2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Asiaticoside B **Chrom-Column:** Arcus C18
Lot Number: 18041322 **Mobile Phase:** ACN: H2O(15%-60%)
Detect Time: 2018-4-13 **Detector:** ELSD
Inject Vol: 20.0 μ L **Tube Temp:** 110 $^{\circ}$ C
Flow: 1.0ml/min **Gas Flow:** 2.9L/min
Sample dissolution MeOH **Impactor:** Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	14.198	392870.938	3123124.500	100.0000