



同田生物
TAUTO BIOTECH

地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0005	Structure 
PRODUCT	Genistein /染料木素	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	18062806	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-06-28	
DATE OF REPORT	2018-06-28	

OTHER NAME 4',5,7-Trihydroxyisoflavone

CHEMICAL FORMULA C₁₅H₁₀O₅

MOLECULAR WEIGHT(MW) 270.2

CAS NUMBER 446-72-0

CHEMICAL FAMILY Flavonoids

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK6—2003	NLT ≥98.0%	99.2%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK6—2003	NA	ND
Water	Q/TFFK6—2003	≤2%	Conforms
Appearance		Pale white crystal	Conforms

STORAGE CONDITIONS

STORAGE Refer to the label

SHELF LIFE 2 Years

Analyst

Checker

Manager

