


PRODUCT ID	E-0187	Structure 
PRODUCT	Ginsenoside-Rc /人参皂昔Rc	
STANDARD TYPE	Qualitative Reference Material/ 含量测定	
LOT NUMBER	18022131	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-02-21	
DATE OF REPORT	2018-02-21	

OTHER NAME	NA
CHEMICAL FORMULA	C ₅₃ H ₉₀ O ₂₂
MOLECULAR WEIGHT(MW)	1079.2685
CAS NUMBER	11021-14-0
CHEMICAL FAMILY	Triterpenes

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK190—2007	NLT ≥98.0%	99.6%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK190—2007	NA	ND
Water	Q/TFFK190—2007	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE Refer to the label

SHELF LIFE 2 Years

Analyst

Checker

Manager

Product Analysis Report


Product Name: Ginsenoside-Rc **Chrom-Column:** C18

Lot Number: 18022131 **Mobile Phase:** ACN:H2O=33:67

Temperature: 25°C **Wavelength:** 203nm

Detect Time: 2018-02-21 **Inject Vol:** 20μL

Analyst: shuzhilin **Flow:** 1ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	8.872	354.090	6767.493	0.3170
2	12.367	103280.188	2128301.750	99.6830