



Mail: sales2@tautobiotech.com export2@tautobiotech.com

CCN1CCC23C4C1CC5=C3C(=C(C=C5)O)C2(C)C=C4


ANALYTICAL RESULT

STORAGE CONDITIONS

Manager

Product Analysis Report

Product Name: Bullatine A **Chrom-Column:** C18
Lot Number: 18022722 **Mobile Phase:** ACN:0.02MNa2HPO4(PH7.3)=6
0:40
Temperature: 25°C **Wavelength:** 210nm
Detect Time: 2018-02-27 **Inject Vol:** 20μL
Sample: MeOH **Flow:** 1ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	4.333	14142.621	155052.938	1.3003
2	4.455	12041.848	89916.383	0.7540
3	5.383	145.721	1728.032	0.0145
4	5.778	9018.822	99254.203	0.8323
5	7.755	412382.156	11525518.000	96.6515
6	8.857	331.522	5463.231	0.0458
7	9.303	211.861	2938.658	0.0246
8	9.835	230.673	3386.977	0.0284
9	10.977	183.723	3906.300	0.0328
10	14.575	347.402	8373.399	0.0702
11	17.678	768.308	18892.301	0.1584
12	18.590	147.147	3656.323	0.0307
13	19.825	132.775	3471.844	0.0291
14	20.673	180.555	3265.530	0.0274