



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0538	Structure 
PRODUCT	Isoschaftoside/异夏佛塔苷	
STANDARD TYPE	Qualitative Reference Material/ 含量测定	
LOT NUMBER	18022831	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-02-28	
DATE OF REPORT	2018-02-28	

OTHER NAME	NA
CHEMICAL FORMULA	C ₂₆ H ₂₈ O ₁₄
MOLECULAR WEIGHT(MW)	564.492
CAS NUMBER	52012-29-0
CHEMICAL FAMILY	Flavonoids

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK573—2008	NLT ≥98.0%	98.6%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK573—2008	NA	ND
Water	Q/TFFK573—2008	≤2%	Conforms
Appearance		Yellow Powder	Conforms

STORAGE CONDITIONS

STORAGE	Refer to the label
SHELF LIFE	2 Years

Analyst

Checker


Manager

Product Analysis Report

Product Name: Isoschaftoside Chrom-Column: C18

Lot Number: 18022831 Mobile Phase: ACN:0.1%H3PO4(15: 85)

Detector: UV Wavelength: 272nm
Detect Time: 2018-02-28 Inject Vol: 20ul
Analyst: shuzhilin Flow: 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	15.048	881.767	23017.500	0.2468
2	16.865	498.447	14201.900	0.1523
3	21.567	238.448	9814.544	0.1053
4	24.302	235011.469	9198239.000	98.6417
5	28.338	1716.727	79630.555	0.8540