


PRODUCT ID	E-0008	Structure 
PRODUCT	Ginkgolide B /银杏内酯B	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	18032622	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-3-26	
DATE OF REPORT	2018-3-26	

OTHER NAME	NA
CHEMICAL FORMULA	C ₂₀ H ₂₄ O ₁₀
MOLECULAR WEIGHT(MW)	424.39856
CAS NUMBER	15291-77-7
CHEMICAL FAMILY	Diterpenoids

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK20—2004	NLT ≥98.0%	100.0%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK20—2004	NA	ND
Water	Q/TFFK20—2004	≤2%	Conforms
Appearance		White Powder	Conforms

STORAGE CONDITIONS

STORAGE	Refer to the label
SHELF LIFE	2 Years

Analyst

Checker

Manager

Product Analysis Report

Product Name: Ginkgolide B

Chrom-Column: Arcus C18

Lot Number: 18032622

Mobile Phase: ACN: H2O(28: 72)

Detect Time: 2018-3-26

Detector: ELSD

Inject Vol: 20.0 μ L


Tube Temp: 110 $^{\circ}$ C

Flow: 1.0ml/min

Gas Flow: 2.5L/min

Sample dissolution MeOH

Impactor: Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	17.228	840292.250	15392192.000	100.0000