



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059

86-21-51320588-8008/8035/8026/8025


传真：86-21-51320502

QQ:800015916

MSN:Tauto_shanghai@hotmail.com

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0010	Structure 
PRODUCT	Ginkgolide J /银杏内酯J	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	18021716	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-02-17	
DATE OF REPORT	2018-02-17	

OTHER NAME	NA
CHEMICAL FORMULA	C ₂₀ H ₂₄ O ₁₀
MOLECULAR WEIGHT(MW)	424.39856
CAS NUMBER	107438-79-9
CHEMICAL FAMILY	Diterpenoids

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK16—2004	NLT ≥98.0%	98.5%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK16—2004	NA	ND
Water	Q/TFFK16—2004	≤2%	Conforms
Appearance		White Crystal	Conforms

STORAGE CONDITIONS

STORAGE	Refer to the label
SHELF LIFE	2 Years

Analyst

Checker

Manager

Product Analysis Report

Product Name: Ginkgolide J

Chrom-Column:

Lot Number: 18021716

Mobile Phase: H₂O: MeOH (70: 30)

Detect Time: 2018-02-17

Detector: ELSD

Inject Vol: 20 µ L


Tube Temp: 110°C

Flow: 1.0ml/min

Gas Flow: 2.9 L/min

Analyst: zhanghong

Impactor: Off



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	10.165	19319.486	217303.813	0.5492
2	10.282	21943.959	344332.844	0.8702
3	15.265	810931.875	39009288.000	98.5807