



地址：上海市张江高科技园区爱迪生路326号 邮编：201203

销售热线：400-820-3059


86-21-51320588-8008/8035/8026/8025

传真：86-21-51320502

QQ:800015916

URL: www.tautobiotech.com

Mail: sales2@tautobiotech.com export2@tautobiotech.com

PRODUCT ID	E-0017	Structure 
PRODUCT	Ginkgolic Acid (C15:1) /银杏酸 (C15:1)	
STANDARD TYPE	Quantitative Reference Material (含量测定)	
LOT NUMBER	18011024	
ASSAY METHOD	HPLC,MS,NMR	
DATE OF SAMPLE	2018-01-10	
DATE OF REPORT	2018-01-10	

OTHER NAME	NA
CHEMICAL FORMULA	C ₂₂ H ₃₄ O ₃
MOLECULAR WEIGHT(MW)	346.5076
CAS NUMBER	22910-60-7
CHEMICAL FAMILY	Phenolic Acid

ANALYTICAL RESULT

TEST	METHOD	SPECIFICATION	RESULT
HPLC	Q/TFFK32—2005	NLT ≥98.0%	98.7%
NMR		Conforms	Conforms
MS		Conforms	Conforms
Residual Solvent	Q/TFFK32—2005	NA	ND
Water	Q/TFFK32—2005	≤2%	Conforms
Appearance		Pale white powder	Conforms

STORAGE CONDITIONS

STORAGE	Refer to the label
SHELF LIFE	2 Years


Analyst

Checker

Manager

Product Analysis Report

Product Name: Ginkgolic Acid (C15:1) **Chrom-Column:** Kromasil C18
Lot Number: 18011024 **Mobile Phase:** MeOH:1%H3PO4 (95:5)
Detector: UV **Wavelength:** 310nm
Detect Time: 2018-01-10 **Inject Vol:** 20ul
Sample dissolution: MeOH **Flow:** 1.0ml/min



No.	Ret.Time min	Height mAU	Area mAU*min	Rel.Area %
1	11.023	3299.081	40868.594	1.2955
2	11.423	169812.156	3113802.500	98.7045